



Dental Equipment Guide

Common Failure • Diagnostic Steps • Quick Fixes

Table of Contents

1.	Imaging & Sensors CBCT, Pano, Intraoral X-Ray, Sensors, Scanner	3
2.	Operatory Equipment Dental Chair, Delivery Unit, Operating Light	4
3.	Sterilization & Infection Control Autoclave, Ultrasonic Cleaner, Washer, Dry Heat	5
4.	Mechanical Systems Compressor, Vacuum Pump, Amalgam Separator	6
5.	Lasers & Surgical Soft/Hard Tissue Lasers, Surgical Drill/Motor	7
6.	Endodontic & Digital Workflow Apex Locator, Rotary Endo Motor, CAD/CAM, Scaler	8
	Quick-Reference Diagnostic Checklist All 22 equipment types at a glance	9
	When to use CapexGPT	10

How to Use This Guide

This guide covers **22 equipment types** across six categories. Each section follows the same structure a field technician uses: **identify symptoms**, **run diagnostics**, and **attempt quick fixes** before scheduling a service call.

When a quick fix isn't enough, that's where CapexGPT comes in. Our Capex Reports give you a data-driven repair vs. replace recommendation with financial analysis — so you never rely on a vendor's word alone.

SAFETY: Always power down and depressurize equipment before performing any diagnostic steps. Follow manufacturer safety protocols and OSHA/CDC guidelines at all times.

1. Imaging & Sensors

Covers: CBCT (Cone Beam 3D), Panoramic X-Ray (2D), Intraoral X-Ray (Wall-Mounted & Handheld), Intraoral Sensors (Digital), Intraoral Scanner

CBCT / Panoramic X-Ray — Unit Won't Complete Scan or Image Is Distorted

Symptoms: Error on startup, arm stops mid-rotation, blurry or streaked 3D/pano images, software won't reconstruct.

- Power cycle the unit completely (full shutdown, wait 60 sec, restart).
- Check for patient positioning errors — movement during scan is the #1 cause of artifacts.
- Verify the chin rest, bite block, and head stabilizers are properly secured.
- Inspect the rotating arm path for physical obstructions or loose covers.
- Check software for error logs — sensor calibration may need to be re-run.
- Confirm the PC meets minimum specs and the connection (USB/Ethernet) is stable.

Intraoral X-Ray (Wall-Mounted & Handheld) — No Exposure or Inconsistent Images

Symptoms: Unit doesn't fire, exposure light/tone absent, image too dark or too light.

- Wall-mounted: Check the circuit breaker on the dedicated outlet and the unit's on/off switch.
- Handheld: Verify battery charge level — low battery = inconsistent output.
- Confirm exposure settings (kV, mA, time) match the anatomy being imaged.
- Test-fire without a sensor to confirm the tube head activates (listen for the tone).
- Check the extension arm for loose joints that could affect beam alignment.

Intraoral Sensors (Digital) — Sensor Not Detected / Image Artifacts

Symptoms: Software shows "no sensor," black or partial images, lines/banding on captures.

- Unplug USB and reconnect to a direct port (avoid USB hubs).
- Check Device Manager (Windows) for driver issues — yellow icon = reinstall driver.
- Test the sensor on a different workstation to isolate PC vs. sensor failure.
- Inspect the cable near the connector — kinks at the junction are the #1 cable failure point.
- Lines or banding = likely internal sensor damage. Visual cracks on the face confirm this.

Intraoral Scanner — Poor Scan Quality or Connectivity Failures

Symptoms: Scans won't stitch, frequent "lost tracking," foggy or incomplete captures, software crashes.

- Clean the scanner tip/mirror with manufacturer-approved wipes (fog = contamination).
- Calibrate the scanner using the provided calibration tool before each session.
- Ensure adequate ambient lighting — too bright or too dark affects optical scanners.
- Update scanner firmware and software to the latest version.
- Check USB-C / proprietary cable for damage; use the original cable only.

PRO TIP: Before assuming a sensor is dead, test on another computer. Roughly 40% of "sensor failure" service calls turn out to be USB port, driver, or software issues.

2. Operatory Equipment

Covers: Dental Chair, Delivery Unit, Dental Operating Light

Dental Chair — Won't Rise, Drifts Down, or No Movement

Symptoms: Chair sinks during procedures, won't lift to full height, motor buzzes but nothing happens.

- Check foot control connection — tubing and wiring secure, unkinked.
- Listen for the hydraulic pump motor. No sound = check the breaker/fuse first.
- Motor runs but no movement = low hydraulic fluid or leaking cylinder seal.
- Inspect the chair base for hydraulic fluid pooling (indicates seal failure).
- Test the foot switch independently — swap with a spare to rule out the control.
- For electric chairs: check the control board for error indicators/LED codes.

PRO TIP: A slowly drifting chair is almost always a worn cylinder seal. Replacement seals typically cost \$150–\$400 vs. \$8,000–\$15,000+ for a new chair.

Delivery Unit — No Water/Air to Handpiece, Weak Flow

Symptoms: Handpiece has no water spray or air, intermittent flow, pressure feels low.

- Check the master water supply valve (under/behind the unit) — confirm it's open.
- Inspect water bottle or direct-connect line for kinks or disconnections.
- Toggle each handpiece holder switch to verify air and water independently.
- Remove the handpiece and test the coupler directly — flow at coupler = handpiece issue.
- Check the inline water filter for sediment clogs (common in older plumbing).
- Verify air pressure at the delivery unit gauge reads 70–80 PSI.

Dental Operating Light — Flickering, Dim, or No Output

Symptoms: Light flickers during procedures, noticeably dimmer than normal, won't turn on.

- Check the power connection at the ceiling arm or wall mount — loose connections cause flickering.
- For LED lights: flickering often indicates a failing driver/power supply board.
- For halogen: replace the bulb first — bulbs dim gradually before full failure.
- Test the on/off switch and dimmer control. Bypass the dimmer if possible to isolate.
- Inspect the articulating arm wiring — repeated repositioning can break internal conductors.
- Check for a reset button on the power module (some models have thermal protection).

SAFETY: Never bypass the check valve on waterlines. Cross-contamination risk violates infection control protocols.

3. Sterilization & Infection Control

Covers: Autoclave (Steam Sterilizer), Ultrasonic Cleaner, Instrument Washer/Disinfector, Dry Heat Sterilizer

Autoclave (Steam Sterilizer) — Won't Reach Temp/Pressure or Wet Packs

Symptoms: Cycle aborts mid-run, error codes (e.g., C699, C501), wet pouches at end of cycle.

- Check water reservoir — refill with distilled water only (tap water = mineral buildup).
- Inspect door gasket for cracks, warping, or debris. Clean with a damp cloth.
- Verify the drain line is clear — disconnect and flush with distilled water.
- Wet packs? Don't overload the chamber — leave space between pouches.
- Check/clean the drain screen and air intake filter.
- Run an empty test cycle to isolate load vs. unit issues.

Ultrasonic Cleaner — Poor Cleaning Performance or No Cavitation

Symptoms: Instruments come out still dirty, no visible cavitation action, unit runs but seems weak.

- Check solution level — the transducer must be fully submerged to create cavitation.
- Replace the cleaning solution (enzymatic solutions lose effectiveness over time).
- Use the foil test: hold aluminum foil in the tank for 30 sec — even pitting = good cavitation.
- Verify temperature setting — most enzymatic solutions work best at 130–150°F.
- Check for solution contamination — heavy debris reduces ultrasonic energy transfer.
- Inspect the tank for cracks or erosion pitting on the bottom surface.

Instrument Washer/Disinfector — Cycle Failures or Incomplete Cleaning

Symptoms: Cycle won't start, instruments still soiled, water not draining, error codes.

- Check water supply and drain connections — kinks or clogs cause most cycle failures.
- Clean the spray arm nozzles — debris blockage is the #1 cause of poor wash performance.
- Verify detergent/rinse agent levels in the dispensing reservoirs.
- Inspect door latch and gasket — an incomplete seal will prevent the cycle from starting.
- Check the drain filter/screen and clean accumulated debris.

Dry Heat Sterilizer — Won't Reach Temperature or Uneven Heating

Symptoms: Takes too long to reach sterilization temp, instruments not sterile per indicators.

- Don't overload — dry heat requires air circulation around every instrument.
- Verify the thermostat setting matches your cycle requirements (typically 320–340°F).
- Check the heating element — no temperature rise at all indicates element failure.
- Confirm the door fully closes and seals (air leaks prevent temp maintenance).
- Use chemical indicators inside the load to verify sterilization, not just external tape.

PRO TIP: Use only distilled water in autoclaves. Tap water mineral deposits are the #1 cause of heating element and reservoir failures — and they're 100% preventable.

4. Mechanical Systems

Covers: Dental Compressor, Vacuum Pump (Suction System), Amalgam Separator

Dental Compressor — Runs Constantly, Won't Build Pressure, or Won't Start

Symptoms: Motor never shutsoff, tankgaugestayslow, audible airleaks, motorhums but won't spin.

Check tank pressure gauge — should cycle between ~90 PSI (cut-in) and ~120 PSI (cut-out).

Listen for air leaks at fittings, regulators, and the tank drain valve.

Drain the tank — moisture accumulation reduces capacity and corrodes internally.

Check the intake filter/muffler for clogs (restricts airflow, reduces output).

Won't start? Check the dedicated circuit breaker and thermal overload reset.

Motor hums but won't spin = possible seized piston or failed start capacitor.

Check the pressure switch — if it doesn't trip at cut-out, it needs recalibration.

PRO TIP: Drain the compressor tank daily. Even oil-free dental compressors accumulate condensation that corrodes the tank and contaminates moisture traps.

Vacuum Pump (Suction System) — Weak Suction or Unusual Noise

Symptoms: HVE barely pulls, saliva ejector has nodraw,pump making screeching/grinding sounds.

Check the vacuum gauge at the pump — should read 5–10 inches Hg.

Inspect and clean chairside trap screens and solids collectors.

Run enzymatic line cleaner through every operator's vacuum lines.

Wet-ring pumps: verify water flow is adequate. Low water = low vacuum.

Dry pumps: check the air filter and impeller housing for debris.

Test suction at each operator individually to isolate one op vs. system-wide.

Screeching = failing bearing or belt. Grinding = debris in impeller. Address immediately.

Amalgam Separator — Full Container, Alarm, or Low Flow

Symptoms: Collectioncanisteralarm, flowrestriction throughthesystem, compliance warning.

Check the collection canister — replace when full per manufacturer schedule.

Inspect inlet/outlet connections for clogs or buildup restricting flow.

Verify the unit is installed in the correct orientation (gravity-fed models are position-sensitive).

Check for air leaks at fittings — leaks reduce separator efficiency.

Confirm your recycling/certification is current — EPA compliance requires proper documentation.

SAFETY: Dental compressors must be on a dedicated circuit. Sharing with other high-draw equipment causes nuisance trips and premature motor failure.

5. Lasers & Surgical Systems

Covers: Laser (Soft Tissue), Laser (Hard Tissue), Surgical Drill/Motor System

Soft Tissue Laser (Diode) — No Output, Weak Cutting, or Fiber Issues

Symptoms: Laser doesn't fire, poor cutting efficiency, fiber tip charred or broken.

- Check the fiber tip — a charred, broken, or improperly initiated tip won't transmit energy.
- Re-initiate (strip and cleave) the fiber per manufacturer instructions.
- Verify power settings match the procedure — too low = poor performance, not a unit failure.
- Inspect the handpiece connector for debris or contamination on the optical interface.
- Check the foot pedal connection — no fire signal = no laser output.
- Confirm the aiming beam is visible — if not, the diode module may need service.

Hard Tissue Laser (Er:YAG / Er,Cr:YSGG) — Reduced Cutting or Water Spray Issues

Symptoms: Laser cutting feels slow, water spray inconsistent, error codes on display.

- Check the water supply — hard tissue lasers require water spray for ablation and cooling.
- Inspect and clean the handpiece tip — debris accumulation reduces beam delivery.
- Verify the air/water lines to the handpiece are connected and free of kinks.
- Check for error codes on the display — most units have built-in diagnostic menus.
- Calibrate the unit per manufacturer schedule — power output drift is common over time.
- Inspect fiber or articulated arm mirrors for contamination (model-dependent).

Surgical Drill / Motor System — Loss of Torque, Speed Fluctuation, or Overheating

Symptoms: Motor stalls under load, speed drops inconsistently, unit feels hot during use.

- Check the handpiece attachment — an improperly seated handpiece reduces torque transfer.
- Verify the irrigation line is flowing — overheating is often lack of coolant, not motor failure.
- Inspect the motor coupling for wear or debris.
- Check speed/torque settings on the control unit — accidental changes are common.
- For implant motors: verify the contra-angle gear ratio matches the programmed setting.
- Lubricate the handpiece per manufacturer specs and test again.

SAFETY: Laser safety protocols are non-negotiable. Ensure all personnel wear proper OD-rated eyewear and that the laser warning sign is posted during any use.

PRO TIP: For diode lasers, 80% of "laser not cutting" complaints are solved by properly re-initiating the fiber tip. Master this step and you'll save most service calls.

6. Endodontic & Digital Workflow

Covers: Apex Locator, Rotary Endodontic Motor, CAD/CAM System, Ultrasonic Scaler

Apex Locator — Erratic Readings or No Signal

Symptoms: Reading jumps erratically, doesn't register the apex, display stays blank.

- Check the lip clip and file clip connections — corroded or loose contacts cause erratic readings.
- Dry the canal adequately — excess sodium hypochlorite or blood gives false readings.
- Replace the battery (often a simple coin cell) — low voltage = unreliable readings.
- Test with a known working file — a bent or fractured file won't read properly.
- Verify the canal has a patent glide path — blocked canals give no reading.
- Clean the clip contacts with alcohol and ensure metal-to-metal contact.

Rotary Endodontic Motor — Auto-Reverse Issues, Torque Errors, or Stalling

Symptoms: Motor reverses unexpectedly, files break more often, unit beeps with torque warnings.

- Verify the torque and speed settings match the file system being used (check the file chart).
- Incorrect settings are the #1 cause of premature file breakage — not motor failure.
- Check the contra-angle head for debris or worn gears (gritty feel when rotated by hand).
- Inspect the file clip/chuck — if it's not gripping properly, the file slips and triggers auto-reverse.
- Test with a new file to rule out a fatigued file vs. motor issue.
- Ensure the handpiece cord isn't being tugged or kinked during use.

CAD/CAM System — Scan Failures, Milling Errors, or Software Crashes

Symptoms: Scan won't stitch, mill produces rough or inaccurate restorations, frequent software errors.

- Scanner: clean the tip, run calibration, and check the USB/connection cable.
- Milling: inspect the burs for wear — dull burs produce rough margins and poor fit.
- Replace milling burs at the manufacturer-recommended interval (don't wait until failure).
- Check the block seating — a loose block in the clamp = inaccurate milling.
- Verify the mill spindle is clean and free of debris buildup.
- Software crashes: check PC RAM and storage — CAD/CAM software is resource-intensive.

Ultrasonic Scaler — Weak Vibration, No Water, or Tip Won't Tighten

Symptoms: Scaling feels weak, water flow absent or inconsistent, insert feels loose.

- Check the power/intensity setting — ensure it's set appropriately for the tip being used.
- Verify the water supply is connected and the flow adjustment is open.
- Inspect the insert tip for wear — replace when 1/2 mm of tip length is lost.
- Magnetostrictive: check the insert stack for corrosion or a cracked ferrite rod.
- Piezo: verify the tip is properly torqued into the handpiece with the manufacturer wrench.
- No vibration at all? Test with a different insert to isolate handpiece vs. insert failure.

PRO TIP: For rotary endo motors, cross-reference your file manufacturer's torque/speed chart every time. Incorrect settings cause more file separations than motor malfunctions.

Quick-Reference Diagnostic Checklist

Print this page and keep it in your equipment room.

Equipment	Check First	Call a Tech If...
CBCT / Pano	Power cycle, patient positioning, software logs	Arm won't rotate, persistent artifacts after recalibration
Intraoral X-Ray	Breaker/battery, exposure settings, test fire	Tube head won't fire after power checks
Intraoral Sensor	USB reconnect, driver check, test on 2nd PC	Cable damage near connector, persistent image lines
Intraoral Scanner	Clean tip, calibrate, update firmware	Hardware error after calibration, broken tip housing
Dental Chair	Foot control, breaker, hydraulic fluid level	Fluid pooling at base, motor runs with no movement
Delivery Unit	Water valve open, check filter, test coupler	No flow after valve/filter check, toggle valve failure
Operating Light	Power connection, bulb/LED, dimmer switch	Flickering after bulb swap (driver board), wiring in arm
Autoclave	Distilled water level, door gasket, drain line	Heating element dead, repeated error codes, BI failures Failed foil test (transducer), cracked tank
Ultrasonic Cleaner	Solution level, temp setting, foil test	Cycle errors after cleaning, pump failure
Washer/Disinfector	Water supply, spray nozzles, door latch	No heat at all (element failure), indicators fail
Dry Heat Sterilizer	Don't overload, thermostat, door seal	Seized piston, tank won't hold pressure (corrosion)
Compressor	Drain tank, intake filter, check for leaks	Grinding from pump, high amp draw, bearing failure
Vacuum Pump	Clean traps, run line cleaner, check gauge	Persistent flow restriction after cleaning
Amalgam Separator	Canister full, inlet/outlet clogs, orientation	No aiming beam, diode module error
Soft Tissue Laser	Re-initiate fiber, check foot pedal, aiming beam	Power output drift, articulated arm damage
Hard Tissue Laser	Water supply, clean tip, check error codes	Motor overheating with coolant flowing, gear wear
Surgical Motor	Handpiece seated, irrigation flowing, settings	Erratic readings with new battery and clean contacts
Apex Locator	Clip contacts, dry canal, replace battery	Motor stalls at correct settings, worn gears
Rotary Endo Motor	Torque/speed settings, check contra-angle, new file	Spindle damage, scanner hardware error, PC failure
CAD/CAM	Calibrate scanner, check burs, block seating	No vibration with multiple inserts, cracked stack
Ultrasonic Scaler	Power setting, water supply, inspect tip wear	

Still Not Sure?

Repair vs. Replace — Answered in 24 Hours.

When a quick fix isn't enough, don't rely on the equipment rep's recommendation. Get a vendor-neutral Capex Report that includes:

CapexScore™ Rating	1-10 score based on age, usage, service history, and reliability.
Financial Analysis	Repair vs. replace cost comparison with Section 179 tax estimates.
Operational Risk	Revenue risk per day of downtime and practice impact assessment.
Lifecycle Benchmarking	Where your equipment stands vs. industry age and parts data.
The Recommendation	A clear repair or replace call. No vendor bias. Just data.

Get Your First Capex Report

CapexGPT.com

(307) 364-2196 | Info@CapexGPT.com

Core \$79.99/mo 1 Report/mo	Pro \$249.99/mo 5 Reports/mo	DSO+ \$599.99/mo 15 Reports/mo
--	---	---

Save 20% with annual plans.

Decide with Data, Not Vendors.